


ORDER FORM		CASE DENTAL & MEDICAL SUPPLY		CONTACT: PAULING TEL: 909.598.0958 / FAX: 909.598.6742			
DR. NAME:			OFFICE NAME:				
ADDRESS:			PHONE:		FAX:		
ITEM	SIZE/COLOR	ITEM #	DESCRIPTION	PKG. QTY	SALE PRICE	QTY	TOTAL
SquirtBrush Refillable							
SquirtBrush Refillable Blister Pack	N/A	RBPC	Refillable Blister Cylindrical + 1 brushtip	1 pc.	\$8.99		
SquirtBrush Refillable Blister Pack	N/A	RBPT	Refillable Blister Taper + 1 brushtip	1 pc.	\$8.99		
SquirtBrush Refillable	0.40 mm	R3SC	Refillable 3S Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.40 mm	R3ST	Refillable 3S Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.48 mm	RSC	Refillable Small Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.48 mm	RST	Refillable Small Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.56 mm	RMC	Refillable Medium Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.56 mm	RMT	Refillable Medium Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.60 mm	RLC	Refillable Large Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
SquirtBrush Refillable Brushtip Blister Pack	N/A	RBPB	Refillable Brushtip x 24	1 pcs.	\$8.99		
SquirtBrush Refillable Brushtip	0.40 mm	RH3SC	Refillable Brush 3S Cylindrical	24/72 pc.	\$12/\$28.8		
	0.40 mm	RH3ST	Refillable Brush 3S Taper	24/72 pc.	\$12/\$28.8		
	0.48 mm	RHSC	Refillable Brush Small Cylindrical	24/72 pc.	\$12/\$28.8		
	0.48 mm	RHST	Refillable Brush Small Taper	24/72 pc.	\$12/\$28.8		
	0.56 mm	RHMC	Refillable Brush Medium Cylindrical	24/72 pc.	\$12/\$28.8		
	0.56 mm	RHMT	Refillable Brush Medium Taper	24/72 pc.	\$12/\$28.8		
	0.60 mm	RHLC	Refillable Brush Large Cylindrical	24/72 pc.	\$12/\$28.8		
SquirtBrush Refillable Brushtip	0.60 mm	RHLT	Refillable Brush Large Taper	24/72 pc.	\$12/\$28.8		
SquirtBrush Professional							
SquirtBrush Professional Straight	0.40 mm	PS3SC	Professional Straight 3S Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.40 mm	PS3ST	Professional Straight 3S Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.48 mm	PSSC	Professional Straight Small Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.48 mm	PSST	Professional Straight Small Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.56 mm	PSMC	Professional Straight Medium Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.56 mm	PSMT	Professional Straight Medium Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.60 mm	PSLC	Professional Straight Large Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
SquirtBrush Professional Straight	0.60 mm	PSLT	Professional Straight Large Taper	24/72/144 pc.	\$54/\$126/\$180		
SquirtBrush Professional Bent	0.40 mm	PB3SC	Professional Bent 3S Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.40 mm	PB3ST	Professional Bent 3S Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.48 mm	PBSC	Professional Bent Small Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.48 mm	PBST	Professional Bent Small Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.56 mm	PBMC	Professional Bent Medium Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
	0.56 mm	PBMT	Professional Bent Medium Taper	24/72/144 pc.	\$54/\$126/\$180		
	0.60 mm	PBLC	Professional Bent Large Cylindrical	24/72/144 pc.	\$54/\$126/\$180		
SquirtBrush Professional Bent	0.60 mm	PBLT	Professional Bent Large Taper	24/72/144 pc.	\$54/\$126/\$180		
SquirtBrush Professional Brushtip Straight	0.40 mm	PHS3SC	Professional Brush 3S Cylindrical	24/72 pc.	\$12/\$28.8		
	0.40 mm	PHS3ST	Professional Brush 3S Taper	24/72 pc.	\$12/\$28.8		
	0.48 mm	PHSSC	Professional Brush Small Cylindrical	24/72 pc.	\$12/\$28.8		
	0.48 mm	PHSST	Professional Brush Small Taper	24/72 pc.	\$12/\$28.8		
	0.56 mm	PHSMC	Professional Brush Medium Cylindrical	24/72 pc.	\$12/\$28.8		
	0.56 mm	PHSMT	Professional Brush Medium Taper	24/72 pc.	\$12/\$28.8		
	0.60 mm	PHSLC	Professional Brush Large Cylindrical	24/72 pc.	\$12/\$28.8		
SquirtBrush Professional Brushtip Straight	0.60 mm	PHSLT	Professional Brush Large Taper	24/72 pc.	\$12/\$28.8		
SquirtBrush Professional Brushtip Bent	0.40 mm	PHB3SC	Professional Brush 3S Cylindrical	24/72 pc.	\$12/\$28.8		
	0.40 mm	PHB3ST	Professional Brush 3S Taper	24/72 pc.	\$12/\$28.8		
	0.48 mm	PHBSC	Professional Brush Small Cylindrical	24/72 pc.	\$12/\$28.8		
	0.48 mm	PHBST	Professional Brush Small Taper	24/72 pc.	\$12/\$28.8		
	0.56 mm	PHBMC	Professional Brush Medium Cylindrical	24/72 pc.	\$12/\$28.8		
	0.56 mm	PHBMT	Professional Brush Medium Taper	24/72 pc.	\$12/\$28.8		
	0.60 mm	PHBLC	Professional Brush Large Cylindrical	24/72 pc.	\$12/\$28.8		
SquirtBrush Professional Brushtip Bent	0.60 mm	PHBLT	Professional Brush Large Taper	24/72 pc.	\$12/\$28.8		
HyssopBrush							
HyssopBrush Blister Pack	0.50 mm	HSB	HyssopBrush Set + 6 Brushtips (3S, S, M)	6/24/144 pc.	\$23/\$60/\$266		
HyssopBrush Brushtip	0.44 mm	HP3SC	HyssopBrush Brush 3S Cylindrical	24/72 pc.	\$7.2/\$21.6		
	0.44 mm	HP3ST	HyssopBrush Brush 3S Taper	24/72 pc.	\$7.2/\$21.6		
	0.50 mm	HPSC	HyssopBrush Brush Small Cylindrical	24/72 pc.	\$7.2/\$21.6		
	0.50 mm	HPST	HyssopBrush Brush Small Taper	24/72 pc.	\$7.2/\$21.6		
	0.56 mm	HPMC	HyssopBrush Brush Medium Cylindrical	24/72 pc.	\$7.2/\$21.6		
	0.56 mm	HPMT	HyssopBrush Brush Medium Taper	24/72 pc.	\$7.2/\$21.6		
	0.60 mm	HPLC	HyssopBrush Brush Large Cylindrical	24/72 pc.	\$7.2/\$21.6		
HyssopBrush Parts	N/A	HPNABC	HyssopBrush Brush Large Taper	24/72 pc.	\$7.2/\$21.6		
HyssopBrush Parts	N/A	HPNABC	HyssopBrush Nozzle Set + Cap	24 pc.	\$48		

SquirtBrush[®]
Solution and Prevention in Oral Health



HyssopBrush[®]
Irrigation Interdental Brush




HyssopBrush
 #HSB

SquirtBrush Refillable




SquirtBrush Refillable
 #RM

SquirtBrush Professional



#PBMC




#PSMC




#RHMC



#RHMT



Bent
 #PHBMC



Straight
 #PHBMC




X-5
 #X5MC


Squirt Brush[®]

Solution and Prevention in Oral Health


Direction:




1. Squeeze and release empty vial in mouthrinse to fill, or use the liquid dropper.




4. Hold bottle in an inverted position. Gently squeeze the bottle while brushing to dispense the content.




2. Insert the base of the brush tip assembly snugly onto the bottle opening.



5. Refill as needed. Rinse bottle and brush after use. Replace brushes weekly or as needed.




3. Gently slide brush between teeth. Never force brush into small/tight spaces.




6. Place cap over brush for storage.

HyssopBrush[®] Irrigation Interdental Brush


Direction:




1. Thread the wire portion of the brush tip through the opening of the brush nozzle to replace the brush tip.




4. Hold bottle in an inverted position. Gently squeeze the bottle while brushing to dispense the content.




2. Secure the brush nozzle by twisting it in a clockwise motion onto the neck assembly. Make sure the tip is positioned entirely into the brush nozzle.



5. Rinse bottle and brush after use. Replace brushes weekly or as needed.



3. Gently slide brush between teeth. Never force brush into small/tight spaces.



6. To refill, simply remove the entire brush from the bottle.

